

DISPOSABLE CRANIAL PERFORATOR

INSTRUCTIONS FOR USE

LANGUAGES

(HU) – Magyar.....	7	(GA) – Gaeilge.....	59
(EN) – English.....	10	(NL) – dutch.....	63
(DE) – Deutsch.....	14	(DA) – Dansk.....	66
(CZ) – Český.....	18	(SV) – svenska.....	70
(SK) – slovenský.....	21	(FI) – Suomalainen.....	73
(RO) – Română.....	24	(EST) – Eesti.....	77
(BG) – Български език.....	28	(LV) – latviešu.....	81
(GR) – Αγγλικά.....	32	(LT) – lietuvių.....	85
(HR) – Hrvatski.....	36	(PL) – Polski.....	88
(SL) – Slovenščina.....	40	(MT) – Inġliż.....	92
(IT) – Italiano.....	43	(RU) – Русский.....	96
(FR) – Français.....	47	(PT) – Português BR.....	100
(ES) – Español.....	51	ی س ر ا ف.....	104
(PT) – Português.....	55	(CH) – 中文.....	108

SZIMBÓLUMOK MAGYARÁZATA | EXPLANATION OF SYMBOLS | SYMBOLER-
KLÄRUNG | VYSVĚTLENÍ SYMBOŮ | VYSVETLENIE SYMBOLOV | EXPLICAȚIA
SIMBOLURILOR | ОБЪЯСНЕНИЕ СИМВОЛОВ | ΕΠΕΞΗΓΗΣΗ ΣΥΜΒΟΛΩΝ | OB-
JAŠNJENJE SIMBOLA | POJASNILO SIMBOLOV | ОБЯСНЕНИЕ НА СИМВОЛИТЕ
| LEGENDA DOS SÍMBOLOS | SPIEGAZIONE DEI SIMBOLI | LEGENDE DES
SYMBOLES | EXPLICACIÓN DE SÍMBOLOS | LEGENDA DOS SÍMBOLOS | MÍNIÚ
AR CHOMHARTHAÍ | UITLEG VAN DE SYMBOLEN | SYMBOLFORKLARINGER
| SYMBOLFÖRTECKNING | SYMBOLIEN SELITYKSET | SÜMBOLITE SELETUS
| SIMBOLU PASKAIDROJUMS | SIMBOLIŲ PAAIŠKINIMAS | OBJAŚNIENIE
SYMBOLI | SPJEGAZZJONI TA' SIMBOLI | ОБЯСНЕНИЕ НА СИМВОЛИТЕ |
ОБЪЯСНЕНИЕ СИМВОЛОВ | LEGENDA DOS SÍMBOLOS | توضیح نمادها |
标识说明

Figyelmeztetés | Warning | Warnung | Varování | Výstraha | Atenție |
Внимание! | Προειδοποίηση | Upozorenje | Opozorilo | Avvertenza |
Avertissement | Advertencia | Aviso | Rabhadh | Waarschuwing | Advarsel | Varning |
Varoitus | Hoiatus | Brīdinājums | Įspėjimas | Ostrzeżenie | Twissija | Внимание! |
Предупреждение | Aviso | اخطار | 提示

Egyszer használatos termék | Disposable product | Einwegprodukt | Likvidace
výrobku | Jednorazový výrobok | Prođus de unică folosință | Изделие за
еднократна употреба | Προϊόν μίας χρήσης | Jednokratni proizvod | Izdelek za
enkratno uporabo | Prodotto utilizzabile una sola volta | Produit à usage unique | Product
desechable | Produto de uso único | Táirge indiúscartha | Wegwerpproduct | Engangsbrug
| Förbrukningsvara | Kertäkäyttöinen tuote | Kõrvaldatav toode | Vienreizlietojamais
produkts | Vienkartinis produktas | Produkt jednorazowego użycia | Prodott disponibbli |
Изделие за еднократна употреба | Продукт одноразового использования | Produto de
uso único | اقلام مصرفی | 一次性产品

Tanulmányozza a használati utasítást! | Consult instructions for use! | Geb-
rauchsanweisung lesen! | Nahlédněte do návodu k použití! | Preštudujte si ná-
vod na použitie! | Consultați instrucțiunile de utilizare! | Прочетете
инструкцията за употреба! | Συμβουλευτείτε τις οδηγίες χρήσης! | Pogledajte upute za
uporabu! | Glejte navodila za uporabo! | Leggere attentamente le istruzioni! | Consulter le
mode d'emploi. | Consultar el manual de instrucciones. | Leia o guia do usuário! | Féach na
teoracha úsáide! | Voor gebruik instructies raadplegen! | Se brugsanvisningen! | Läs
igenom användarinstruktionerna! | Katso käyttöohjeet! | Tutvuge kasutusjuhendiga! |
Skatīt lietošanas instrukcijas! | Vadovaukitės naudojimo instrukcijomis | Przeczytaj
instrukcję użycia! | Ikkonsulta struzzjonijiet għall-użu! | Прочетете инструкцията за
употреба! | Изучите инструкцию по эксплуатации! | Leia o guia do usuário! |
دستور العمل جهت استفاده | 请仔细阅读使用说明书!

STERILE R

Sterilizált termék, sterilizálási eljárás: gamma besugárzás | Sterilized using irradiation | Durch Bestrahlung sterilisiert | Sterilizováno pomocí iradiace | Sterilizované ožarovaním | Sterilizat cu radiații | Стерилизовано изделие, метод на стерилизиране: облъчване с гама-лъчи | Αποστειρώθηκε με ακτινοβόληση | Sterilizirano zračenjem | Sterilizirano z obsevanjem | Prodotto sterilizzato, processo di sterilizzazione utilizzato: raggi gamma | Produit stérilisé, processus de stérilisation: irradiation aux rayons gamma | Producto esterilizado, procedimiento de esterilización: irradiación gamma | Produto esterilizado, proceso de esterilização: irradiação gama | Steirilthe le hionradaíocht | Gesteriliseerd door straling | Steriliseret med stråling | Har strålsteriliserats | Steriloidaan säteilyttämällä | Steriliseeritud, kasutades kiiritust | Sterilizēts, izmantojot apstarošanu | Sterilizuojama naudojantis spinduliuote | Sterylizowane promieniowaniem | Steralizat bl-
užu tal-irradjazzjoni | Стерилизирано изделие, метод на стерилизиране: облъчване с гама-лъчи | Стерилизованный продукт, стерильный процесс: гамма-облучение | Produto esterilizado, proceso de esterilização: irradiação gama | استريل توسط تابش اشعه | 消毒过的产品, 消毒过程: 辐射

Nem újratesterilizálható termék | Do not resterilize! | Nicht erneut sterilisieren! | Nesterilizujte! | Opakovane nesterilizujte! | Nu resterilizați! | Изделието не може да бъде стерилизирано повторно | Μην επαναασταθεροποιείτε! | Ne sterilizirati! | Ne sterilizirajte! | Prodotto non risterrilizzabile | Produit non restérilisable | No reesterilizar | Produto não reesterilizável | Ná hathsterilítear! | Niet nogmaals steriliseren! | Steriliser ikke igen! | Steriliser inte igen! | Älä steriloι uudelleen! | Mitte uuesti steriliseerida! | Nesterilizēt atkārtoti! | Nesterilizuokite! | Nie sterylizować drugi raz! | Tisteralizax mill-ǵdid! | Изделието не може да бъде стерилизирано повторно | Продукт, не подлежащ на повторна стерилизация | Produto não reesterilizável | استريل مجدد نکنید | 不可重新消毒的产品

Sérült vagy nyitott csomagolás esetén a termék nem használható! | Do not use if package is damaged! | Nicht verwenden, wenn Packung beschädigt ist! | Výrobek nepoužívejte, pokud je obal poškozen! | Nepoužívajte, ak je obal poškodený! | A nu se utiliza dacă ambalajul este deteriorat! | Изделието не може да бъде използвано в случай на повредена или отворена опаковка! | Μην το χρησιμοποιείτε εάν το πακέτο είναι κατεστραμμένο! | Ne koristiti ako je paket oštećen! | Ne uporabljajte, če je paket poškodovan! | In caso di imballaggio danneggiato od aperto, il prodotto non è utilizzabile. | Ne pas utiliser en cas d'emballage endommagé ou ouvert. | En caso de envase abierto o dañado el producto no debe utilizarse. | O produto não pode ser usado se a embalagem estiver danificada ou aberta! | Ná húsáidtear má tá damáiste déanta don phacáiste! | Niet gebruiken bij beschadigde verpakking! | Må ikke anvendes, hvis emballagen er beskadiget! | Använd inte om förpackningen är skadad! | Älä käyttää, jos pakkaus on vaurioitunut! | Mitte kasutada, kui pakend on kahjustunud! | Nelietot, ja iepakojumš ir bojāts! | Nenaudokite, jei pažeista pakuotė | Nie używać w razie uszkodzenia opakowania! | Tusazx jekk il-pakkett huwa | Изделието не може да бъде използвано в случай на повредена или отворена опаковка! | Продукт нельзя использовать, если упаковка повреждена или открыта! | O produto não pode ser usado se a embalagem estiver danificada ou aberta! | در صورت آسیب دیدگی جعبه از آن استفاده نکنید

Tárolási hőmérséklet | temperature limit | Temperaturgrenze | Teplotní limit | Teplotný limit | Limita de temperatură | температурен предел | Όριο θερμοκρασίας | Ograničenje temperature | Omejitev temperature | limite di temperatura | limite de température | limite de temperatura | limite de temperatura | Teorainn teochta | Temperatuurlimiet | Temperaturbegrænsninger | Temperaturbegrænsning | Lämpötilarajoitus | Temperatuuripiirang | Temperatūras robeža | Temperatūros riba | Ograniczenie temperatury | Limitu tat-temperatura | температурен предел | Температура хранения | limite de temperatura | حد دما | 温度极限

Relatív páratartalom | humidity limitation | Feuchtegrenzwert | Omezení vlhkosti | Obmedzenie vlhkosti | Limita de umiditate | ограничаване на влажността | Περιορισμός υγρασίας | Ograničenje vlažnosti | Omejitev vlažnosti | limitazione dell'umidità | limitation d'humidité | limitación de humedad | limitação de umidade | Teorainn maidir le bogthaise | Vochtigheidslimiet | Fugtighedsgrænse | Fuktighetsbegrænsning | Kosteuden rajoitus | Niiskuspiirang | Mitruma ierobežojums | Drėgmės apribojimas | Ograniczenie wilgotności | Limitazzjoni tal-umdità | ограничаване на влажността | Относительная влажность | limitação de umidade | محدودیت رطوبت | 湿度限制

93/42/EEC direktiva szerinti meg-felelősség | CE Certification | CE-Zertifizierung | CE Certifikace | CE certifikácia | Certificare CE | съответствие според директива 93/42/ЕИО | Πιστοποίηση CE | CE certifikat | Certifikat CE | Conforme all direttiva 93/42/EEC | Conforme à la directive 93/42/CEE | Conformidad con la directiva 93/42/EEC | Conformidade com a diretiva 93/42/CEE | Deimhniúchán CE | CE-certificatie | CE certifikat | CE Certifiering | CE-sertifiointi | CE-sertifitseerimine | EK sertifikācija | CE atestavimas | Certyfikat CE | Ćertifikazzjoni CE | съответствие според директива 93/42/ЕИО | Соответствие Директиве EEC/93/42 | Conformidade com a diretiva 93/42/CEE | گواهی CE 93/42/EEC | 符合93/42/EEC指示

Lejárati dátum: elérése esetén a termék nem felhasználható! | Use-by date | Haltbarkeitsdatum | Použitelnost do | Dátum použiteľnosti | Data limită de utilizare | Срок на годност: след изтичането му изделието не може да се използва! | Χρήση με βάση Ημερομηνία | Rok upotrebe | Datum uporabe | Data di scadenza: in seguito al raggiungimento della data di scadenza il prodotto non è più utilizzabile! | Date d'expiration: si elle est atteinte, le produit ne peut plus être employé! | Fecha de caducidad: no usar el producto después de la fecha de caduciad. | Data de validade: o produto não pode ser usado após o vencimento! | Dáta éagtha | Gebruiken voor datum | Anvendes inden | Sista förbrukningsdag | Käytä päivämäärään mennessä | Aegumiskuupäev | Lietošanas termiņa beigu datums | Naudokite iki nurodyto termino | Data ważności | Užu sad-data | Срок на годност: след изтичането му изделието не може да се използва! | Срок годности: по истечении срока годности продукт использовать запрещено! | Data de validade: o produto não pode ser usado após o vencimento! | تاريخ مصرف | 保质期过期或不许使用产品!

REF

Termékazonosító kód | Catalogue number | Katalognummer | Katalogové číslo | Katalógové číslo | Număr de catalog | Идентификационный код | Αριθμός καταλόγου | Kataloški broj | Kataloška številka | Codice di identificazione del prodotto | Code d'identification du produit | Número de identificación del producto | Código de identificação do produto | Uimhir chatalóige | Catalogusnummer | Bestellingsnummer | Katalognummer | Luettelonumero | Kataloginumber | Kataloga numurs | Katalogo numeris | Numer katalogowy | Numru tal-katalogu | Идентификационный код | Идентификационный код продукта | Código de identificação do produto | شماره کاتالوگ | 产品识别码

LOT

Termékszéria jelölése | Batch code | Batch-Code | Číslo dávky | Kód šarže | Codul șarjei | Серия на изделието | Κωδικός παρτίδας | Šifra serije | Koda serije | Indicazione della serie di produzione del prodotto | Marquage série produit | Código de lote | Marcação da série do produto | Cód baise | Batchcode | Batchkode | Satsnummer | Eränumero | Partiitähis | Partijas kods | Partijos kodas | Kod partii | Kodići tal-bačić | Серия на изделието | Маркировка серии продуктов | Marcação da série do produto | كد محصول | 产品型号

A gyártó cég megnevezése | Manufacturer | Hersteller | Výrobce | Výrobca | Producător | Название компании-производителя | Κατασκευαστής | Proizvođač | Proizvajalec | Denominazione della società produttrice | Désignation du fabricant | Empresa fabricante | Nome da empresa | Monaróir | Fabrikant | Fabrikant | Tillverkare | Valmistaja | Tootja | Ražotājs | Gamintojas | Producent | Manifattur | Название компании-производителя | Производител | Nome da empresa fabricante | سازنده | 制造企业名称

Tartsa szárazon | Keep dry | Trocken halten | Uchovávejte v suchu | Uchovávejte v suchu | A se păstra uscat | Пази сухо | Διατηρήστε το στεγνό | Čuvati na suhom | Hraniti na suhem | Mantenere asciutto | Garder au sec | Mantener seco | Manter seco | Coinnigh an fheiste tirim | Drooghouden | Opbevaars tørt | Förvaras torrt | Pidā kuivana | Hoida kuivana | Uzglabāt sausā vietā | Laikyti sausoje vietoje | Przechowywać w suchym miejscu | Žomm xott | Держите в сухости | Manter seco | خشک نگه دارید | 保持干燥

Napfénytől óvja | Keep away from sunlight | Vor Sonnenlicht schützen | Chraňte před slunečním zářením | Uchovávejte mimo dosahu slnečného svetla | A se păstra departe de lumina soarelui | Да се пази от слънчева светлина | Φυλάσσετε μακριά από το φως του ήλιου | Čuvati dalje od sunčeve svjetlosti | Hraniti na temnem mestu | Tenere lontano dalla luce solare | Tenir à l'écart de la lumière du soleil | Mantener alejado de la luz solar | Afaste-se da luz solar | Coinnigh amach ó sholas na gréine | Uit de buurt van zonlicht houden | Opbevaars væk fra sollyls | Förvaras skyddad från solljus | Suojaa auringonvalolta | Hoida eemal päikesevalgusest | Sargāt no saules gaismas | Laikyti sausoje vietoje | Przechowywać z dala od światła słonecznego | Žomm 'il bogħod mid-dawl tax-xemx | Беречь от солнечного света | Afaste-se da luz solar | از نور خورشید دور باشید | 远离阳光

Rx Only

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NLO-11/14-3.0, NLO-9/13-3.0, NLO-7/11-3.0, NLO-6/9-3.0

NLO-11/14-1.5, NLO-9/13-1.5, NLO-7/11-1.5, NLO-6/9-1.55

**Important information,
please read before use!**

**Failure to read and follow instructions for
this device may result in serious patient injury or death.**

**Single use product! It is intended for use on
a single patient during a single procedure.**

ABBREVIATION:

CDS Process of Cleaning – Disinfection – Sterilization

DELIVERY CONTENT:

The articles NLO-x/y-z NeuroLine single use cranial perforators are individually wrapped and delivered in sterile condition, ready for use.

**Use only originally wrapped,
sterile perforators!**

The perforator must be considered non-sterile and therefore it is prohibited to use if

- the Gamma indicator symbol on the outer packaging is not red or missing, or
- the expiry date has passed, or
- the blister pack is damaged or open.

Repeated use of disposable product carries a risk of infection. Never use the perforator on more than one patient!

The procedure CDS damages the perforator, which seriously endangers the patient in case of repeated applications, as well as it terminates the manufacturer's liability with immediate effect.

OPERATING CONDITIONS:

- Temperature limit: +10 to +40 °C
- Humidity limitation: 10 – 90 %

TRANSPORT AND STORAGE CONDITIONS:

- Temperature limit: -10 to +50 °C
- Humidity limitation: 10 – 90 %

OPERATING CONDITIONS:

1. Pairing with universal Hudson fitting. Make sure that the perforator shaft is completely and safely implanted into the drive.
2. Use the drive in a clockwise direction and „forward“ mode.
3. The recommended speed of use is 800-1000rpm (up to 1300rpm). Use a drive only of the above speed range.

INTENDED USE:

NeuroLine perforators are designed to make holes in the skull bone the way that the perforator's integrated safety mechanism releases automatically after passing through the bone and stops drilling.

INTENDED USER:

NeuroLine perforators can only be used by neurosurgeons who have the required skill and experience.

FIELD OF APPLICATION:

The NLO-x/y-3.0 models are for use on skulls at least 3.0mm thick.

The NLO-x/y-1.5 models are for use on thin skulls, respectively on thin areas of the skull (e.g. area of the temple or nape) if the skull is at least 1.5mm thick.

It is the user's responsibility to act with adequate care when choosing the suitable type of perforator, considering the patient and its health condition and also the operation areas of the skull. Any selection has been made for that patient, careful use is essential when using the chosen perforator.

Function examination

Check the perforator before surgical application!

The perforator cutting edges are very sharp; therefore, protect your fingers with a piece of multilayer gauze while performing the tests.

Steps of checking operability:

1. Hold the perforator as shown in the figure 1. and press the drill-bit. Expectation: smooth, firm springing effect.
 2. Hold the perforator as shown in the figure 2. and – while grabbing the drive shaft – rotate the drill parts (inner and outer drills). Expectation: unobstructed rotation.
- If any of the checking steps do not result in the indicated test output, stop using the perforator, return it to the manufacturer in its original packing. Do not apply any CDS procedures.

SURGICAL APPLICATION:

1. To start perforation, place the non-skid tip of the perforator on the surface of the skull, hold perpendicular and press firmly onto the bone, then drill with the appropriate speed
2. If the compressive force is terminated and/ or the perforator is not driven during the drilling period, then the drill stops. If necessary, the drilling can be restarted at any time. If you notice any abnormality, do not continue to use the perforator.
3. When the perforator reaches the inner surface of the skull, loss of pressure occurs at the perforator drill head, then the perforator safety mechanism is activated and the drill is stopped: the hole is completed.
4. NeuroLine perforators form a bone plate at the bottom of the hole (figure 3.), which protects the dura from the cutting edges and shifts the dura from the bone. The bone plate can be removed simply by using a forceps or an elevator. If you suspect that the dura adheres to the bone (e.g. according to the cranial suture or in cases of older patients), be especially careful to avoid damage of the dura.

Warnings

1. The surgeon should always be aware that it may happen that the perforator's safety mechanism does not release. Avoid excessive pressing force to prevent patient injury, the damage of the dura and penetration into the brain.
2. Reduce the pressing force and proceed with caution when the perforator is approaching the inner surface of the skull, dura or brain injury may occur.
3. Hold the driver unit and the perforator perpendicular (90°) to the skull. The perforator is optimized for stable directional drilling angle of 90° degrees (figure 4.). The perforator's design is optimized for stable directional drilling angle of 90° degrees, do not rock or tilt it in order to speed up the drilling process. Failure to observe this instruction could cause product malfunction, patient injury or dural laceration.
4. Extreme caution must be taken when perforating in areas of the skull that have 1mm or greater variations in bone thickness. Drill may nick or cut dura or brain (similar effect as not being 90°).
5. A condition such as the dura adheres, or high intracranial pressure or any other underlying disorders exist at the area of penetration, the perforator may nick or cut the dura or the brain.
6. Take special care NOT to exceed the speed of 1300rpm. This may result in damage of patient's dura or brain.
7. The NeuroLine perforators are designed to release safely and reliably, but caution must be observed at all times, in particular at the areas of thin-boned, for example the area

of the temporal bone, or for children, newborns, elderly or patients with sick/damaged bone. Drilling should be proceeded with extreme caution, because bone consistencies, strength and thickness may be different in these cases, and may result in injury to the patient: damage of the dura or the brain may occur.

8. Not recommended for use the perforator in a previously drilled burr hole or near to such areas, because the risk of injury is significantly higher: damage of the dura or the brain may occur.
9. Perforators can only be used with motor systems designed for medical use that have Hudson-end, because the use of any other systems may result in damaging the patient's dura or brain.

LIMITED WARRANTY:

The Manufacturer grants that this product is free from defects regarding materials and workmanship, and – if the conditions and circumstances described in this guide are fulfilled –, safely and optionally functioning for use on a single patient during a single procedure.

If this product fails to meet this warranty, the Manufacturer will, at its discretion, replace it or refund purchase price.

In case any problem occurs during the operation relating to the device, it must be returned to the manufacturer for investigation. Autoclave sterilization makes the investigation process much harder, or even impossible. Therefore, we expressly recommend making disinfection wash only. The device must be put in a proper protecting package to avoid risk of infection during transport. Regarding the return process, please contact the manufacturer or distributor.

This warranty is void, if

- the product is not used in accordance with its intended purpose and/or as described in this instruction guide, or
- damage occurred to the product, or
- the customer failure notification announced after the expiry date, or
- used any phase of process CDS, including gamma radiation, or
- the product is repeatedly used in different surgical procedure.

Customer complaint by the manufacturer is to be investigated only if none of the factors causing the loss of validity of the warranty factor exists, as well as the original manufacturer's packaging remains at disposal.

Any other expressed or presumptive warranties are hereby disclaimed.

DISPOSAL OF PRODUCT:

After use, the institutional and/or in accordance with national legislation and standards for the treatment of infected devices taking into account the waste treated.

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